

Fume Hood Alarm, Issue, and Best Management Practice Responsibilities

Call IHS for initial response to ANY fume hood issue at (706) 721-2663

Lab Staff are responsible for the following annual certification correction items:

- 1. Glass sash should be kept closed when not in use.
- 2. Loose paper towels and KimWipes should not be left in hood.
- 3. Low rear baffle opening should not be blocked.
- 4. Chemical residue should not be on hood surfaces.
- 5. Hood storage should be minimal.
- 6. Large equipment should be raised 1" to 2" (inches) off base.
- 7. All items should be inset 6" (inches) from glass sash opening.

Lab Equipment Services (LES) is responsible for the following annual certification correction items: Contact at (706) 721-6124 or (706) 721-3274

- 8. Baffle adjustment should work properly.
- 9. Sash should move easily and not drift.
- 10. Low airflow alarm should work properly (properly calibrated).
- 11. Lights should work adequately.

Facilities (Work Management) is responsible for the following annual certification correction items: Contact at www.www.www.ement) or (706) 721-2434

- 12. Hood outlet should be properly wired and grounded.
- 13. Noise levels should be acceptable.
- 14. Interior side panels should be secured.

Additional Information

- Work management is responsible for the function/operation of the fume hood.
- Phoenix Alarm Systems fumes hoods are calibrated by outside contractor, submit work order to facilities.